



SRMS Directory

Accreditation number: SRMS 0001

International standard: ISO 17034:2016
Swiss standard: SN EN ISO 17034:2017

Sigma-Aldrich Production GmbH
Industriestrasse 25
9471 Buchs

Head: Michael Mühlhofer
Responsible for MS: Dr. Andrea Schwaiberger
Telephone: +41 81 755 2852
E-Mail: andrea.schwaiberger@merck-group.com
Internet: <http://www.sigmaaldrich.com>
Initial accreditation: 14.11.2007
Current accreditation: 05.09.2022 to 04.09.2027
Scope of accreditation see: www.sas.admin.ch
(Accredited bodies)

Scope of accreditation as of 05.01.2024

Producer of certified reference materials in the fields of organic compounds, organic and inorganic pure substances, solutions, microbiological discs and matrix reference materials

Reference material matrix / artefact	Properties characterized	Type	Approach used to assign property values
CHEMICAL REFERENCE MATERIALS			
Organic compounds	Mass fraction of main or minor component	CRM	Quantitative determination by quantitative performance NMR spectroscopy (qNMR) (a, b, d) ¹
Organic compounds	Mass fraction of main or minor component	CRM	Quantitative determination by LC-MS (IDMS) (a, d) ¹
Organic compounds	Mass fraction of main or minor component	CRM	Quantitative determination by LC-MS (a, d) ¹
Organic compounds	Mass fraction of main or minor component	CRM	Quantitative determination by LC-CAD (a, d) ¹
Organic compounds	Mass fraction of main or minor component	CRM	Quantitative determination by LC-UV (a, d) ¹
Organic compounds	Mass fraction of main or minor component	CRM	Quantitative determination by GC-MS (IDMS) (a, d) ¹



SRMS Directory

Accreditation number: SRMS 0001

Reference material matrix / artefact	Properties characterized	Type	Approach used to assign property values
Organic compounds	Mass fraction of main or minor component	CRM	Quantitative determination by GC-MS (a, d) ¹
Organic compounds	Mass fraction of main or minor component	CRM	Quantitative determination by GC-FID (a, d) ¹
Organic and inorganic pure substances	Mass fraction of the main component	CRM	Quantitative determination by metrological precision titration (a, d) ¹
Organic and inorganic pure substances	Mass fraction of main analyte	CRM	Quantitative determination by ICP-OES (a, d) ¹
Organic and inorganic pure substances	Mass fraction of minor components in the pure substance	CRM	Quantitative determination by ICP-OES, ICP-MS, AAS (a, d) ¹
Organic and inorganic pure substances	Mass fraction of minor components in the pure substance	CRM	Quantitative determination by IC (a, d) ¹
Organic and inorganic pure substances	Mass fraction of main analyte	CRM	Quantitative determination by IC (a, d) ¹
Solutions	Mass fraction of single or multiple components	CRM	Gravimetric production with high precision weighing, starting from pure or high purity starting materials, the content of which has been characterized by the 100% minus impurities approach or by qNMR, LC-MS, LC-CAD, LC-UV, LC-IDMS, IC, titration, GC-MS, GC-IDMS, GC-FID or ICP-OES. (b, e) ¹
Solutions	Mass fraction of single component	CRM	Quantitative determination of the main analyte by metrological precision titration (a, d) ¹
Solutions	Mass fraction of single component	CRM	Quantitative determination of the main analyte by ICP-OES (a, d) ¹
Solutions	Mass fraction of single or multiple components	CRM	Quantitative determination of minor components by ICP-OES, ICP-MS, AAS (a, d) ¹
Solutions	Mass fraction of single or multiple components	CRM	Quantitative determination of the main analyte and/or minor components by IC (a, d) ¹
Solutions	Mass fraction of single or multiple components	CRM	Quantitative determination of the main analyte and/or minor components by quantitative NMR spectrometry (qNMR) (a, d) ¹
Solutions	Mass fraction of single or multiple components	CRM	Quantitative determination of the main analyte and/or minor components by LC-MS (IDMS) (a, d) ¹



SRMS Directory

Accreditation number: SRMS 0001

Reference material matrix / artefact	Properties characterized	Type	Approach used to assign property values
Solutions	Mass fraction of single or multiple components	CRM	Quantitative determination of the main analyte and/or minor components by LC-MS (a, d) ¹
Solutions	Mass fraction of single or multiple components	CRM	Quantitative determination of the main analyte and/or minor components by LC-CAD (a, d) ¹
Solutions	Mass fraction of single or multiple components	CRM	Quantitative determination of the main analyte and/or minor components by LC-UV (a, d) ¹
Solutions	Mass fraction of single or multiple components	CRM	Quantitative determination of the main analyte and/or minor components by GC-MS (IDMS) (a, d) ¹
Solutions	Mass fraction of single or multiple components	CRM	Quantitative determination of the main analyte and/or minor components by GC-MS (a, d) ¹
Solutions	Mass fraction of single or multiple components	CRM	Quantitative determination of the main analyte and/or minor components by GC-FID (a, d) ¹
MICROBIOLOGICAL REFERENCE MATERIALS			
Plano-convex discs	Number of colony-forming bacteria, yeasts, moulds	CRM	Cultural quantitative methods of determination by a laboratory using reference methods (a) ¹
FOOD REFERENCE MATERIALS			
Milk powder	Mass fraction of main or minor components	CRM	Quantitative determination by LC-MS (IDMS) (a, d) ¹

In case of contradictions in the language versions of the directories, the German version shall apply.

¹ Categories according to the standard SN EN ISO 17034:2017 § 7.12.3:

- using a single reference measurement procedure (as defined in ISO/IEC Guide 99) in a single laboratory
- characterization of a non-operationally defined measurand using two or more methods of demonstrable accuracy in one or more competent laboratories
- characterization of an operationally defined measurand using a network of competent laboratories
- value transfer from an RM to a closely matched candidate RM performed using a single measurement procedure performed by one laboratory
- characterization based on mass or volume of ingredients used in the preparation of the RM



SRMS Directory

Accreditation number: SRMS 0001

Abbreviation	Signification
CAD	Charged aerosol detector
CFU	Colony forming unit
CRM	Certified reference material
GC	Gas chromatography
IC	Ion chromatography
ICP-OES	Inductively coupled plasma optical emission spectrometry
ICP-MS	Inductively coupled plasma combined with mass spectrometry
IDMS	Isotope dilution mass spectrometry
LC	Liquid chromatography
MS	Mass spectrometry
NMR	Nuclear magnetic resonance
UV	Ultraviolet

* / * / * / * / *